





Omeprazole Related Substances (USP)

On Endurus[®] C8 Classic 5 μ m 250X4.6



Gradient Table

Time (Min)	MPA	MPB
0	75	25
12	75	25
22	50	50
45	50	50
45.1	75	25
55	75	25

Test Condition (As USP)

Column:

Endurus[®] C8 Classic 250x4.6 mm, 5 μ m

Injection: 40 μ L

Detection: UV 264 nm

Flow Rate: 0.8 mL/min

Mobile Phase:

B: Acetonitrile

A: Dissolve 0.725 g of monobasic sodium phosphate and 4.472 g of anhydrous dibasic sodium phosphate in 1000 mL of water. Dilute 250 mL of this solution with water to 1000 mL. If necessary, adjust with phosphoric acid to a pH of 7.0.

Diluent: Acetonitrile and Buffer (1:3)

Temperature: Ambient

Autosampler temperature: 5 $^{\circ}$ C

Gradient: See Table

Sample preparation:

0.6 mg/mL of Omeprazole in Diluent.

Chromatographic data

No.	Compound	Retention Time (min)	Tailing Factor	RRT
1	Desmethoxy omeprazole	21.69	0.99	0.95
2	Omeprazole	22.78	1.02	1
3	Omeprazole thioxo	35.655	1.02	1.57

omeprazole exerts its therapeutic effects by selectively and irreversibly inhibiting the H⁺/K⁺-ATPase system in gastric parietal cells. This inhibition leads to a reduction in stomach acid secretion, alleviating symptoms associated with acid-related disorders. Understanding the biochemical mechanism of omeprazole's action provides insights into its clinical utility and underscores its role as a cornerstone therapy for acid-peptic diseases.

As per United States Pharmacopeia

Suitability requirements

Tailing factor: NMT 1.5

Relative retention time (RRT): 1.0 for omeprazole

ISO 9001
CERTIFIED

To place an order, call +91-90811 21133, or contact your local Force Scientific distributor:

Force Scientific • 207 & 210-213 The Space • Kunal Cross Road • Gotri, Vadodara, GJ 390021 • www.forcescientific.com

Endurus[®] is a registered trademark of Force Scientific.